

Tape plaster RiverPlast



Technical characteristics:

Classic type

- Made on perfect woven 100% cotton base.
- Optimum adhesiveness.
- Can be removed painlessly from skin. Don't leave glue traces on skin.
- Hypoallergenic.
- Microporous structure permits skin to "breath".

Transparent type

- Waterproof. Made on polyethylene base.
- Indispensable in case of contact with water.
- Can be easily torn off in any direction without usage of scissors.
- Invisible on skin and hypoallergenic.
- Microporous structure permits skin to "breath".

Ultrapore type

- Made on non-woven base.
- Hypoallergenic, recommended for sensitive and child's skin.
- Can be used instead of wound plaster in case of minor cuts, abrasions and scratches.
- Almost invisible due to skin color of the plaster base.
- Can be easily torn off in any direction without usage of scissors.
- Microporous structure permits skin to "breath".

Sizes of plaster: 1x500 cm, 2x500 cm, 3x500 cm, 5x500 cm.
Available in handy plastick or ecomony paper packing.

Usage:

For additional fixation of bandages, catheters, infusion/transfusion sets and other medical devices.

Marking and labeling are made in Ukrainian and comply with all requirements of standards.

The product is registered in Ukraine and allowed to use in medical practice.

Wound plaster «IGAR»



Technical characteristics:

Classic type

- Made on natural 100% cotton base.
- Hypoallergenic.
- Are recommended for people with sensitive skin.
- Microporous structure permits skin to "breath".
- Almost invisible due to skin color of the plaster base.

Waterproof type

- Made on medical PVC base.
- Waterproof. Good adhesiveness even when wet.
- Low allergenic adhesive film doesn't leave traces on skin.
- Microporous structure permits skin to "breath".
- Almost invisible due to skin color of the plaster base.

Elastic type

- Made on waterproof elastic base.
- Hypoallergenic.
- Perfect conformity due to elasticity.
- Ideal for joint areas.
- Microporous structure permits skin to "breath".
- Almost invisible due to skin color of the plaster base.

Sizes of wound plasters: 7,2x1,9cm и 7,6x2,5 cm.
Sterilized with ethilene oxide gas.

Usage:

First aid in case of minor cuts, abrasions and scratches, on injection sites.

Marking and labeling are made in Ukrainian and comply with all requirements of standards.

The product is registered in Ukraine and allowed to use in medical practice.

**A doctor should fix
his mind on curing**



IGAR Ltd. www.igar.ua
20, Metalistiv str., Kiev, 03057
phone/fax: 044-585-34-55