


## Foley Catheter «IGAR»

For enduring drainage  
of urinary bladder



### Technical Characteristics:

- Made of latex, with silicon coating, which provides good compatibility with human tissue, protects from irritations and allergic reactions and eases passage through urethra.
- The symmetrical balloon permits to place catheter tip accurately and safely.
- Two side apertures on the closed end of catheter.
- Balloon of 30 ml.
- Rounded form and reinforced connective parts provide optimal strength.
- Air valve (check valve) permits all commonly used syringes for easy inflation and deflation.
- The thick and firm catheter shaft provides safe connection with drainage pipe.
- Catheter is optimally rigid; it becomes softer at body temperature what reduces possible discomfort.
- Apertures in proper proportion are designed for effective drainage.
- Standard color coding.
- Sterilized with ethylene oxide gas.

### Usage:

- Drainage on continuing basis of urinary bladder in case of difficult or impossible urination, unconscious state or in early postoperation period.
- Urinary bladder irrigation.
- Kidney function control.
- Incontinence of urine.

Color code	Gauge/Size
Green	14 (4,7 x 390mm)
Orange	16 (5,3 x 390mm)
Red	18 (6,0 x 390mm)
Yellow	20 (6,7 x 390mm)
Violet	22 (7,3 x 390mm)
Blue	24 (8,0 x 390mm)
Rose	26 (8,7 x 390mm)

Marking and labeling are made in Ukrainian and comply with all requirements of standards.

The product is registered in Ukraine and allowed to use in medical practice.

## Nelaton Catheter «IGAR»

For short-term catheterization  
of urinary bladder



### Technical Characteristics:

- Made of transparent non-toxic medical PVC for minimal irritations and allergic reactions during catheterization.
- Thermoplastic catheter material becomes softer at body temperature what reduces possible discomfort.
- No need in additional lubricants due to smooth and silky surface.
- Two side apertures on the distal end of catheter for effective drainage without risk of obturation.
- Distal end of catheter is rounded and perfectly finished to provide non-traumatic catheterization.
- Connector allows usage with urinary bags of any type.
- Length of catheter: female - 18±2 см; male - 40±2 см.
- Standard color coding.
- Sterilized with ethylene oxide gas.

### Usage:

- Short-term catheterization of urinary bladder
- Instillation of medical solutions into urethra and/or urinary bladder.

Цветовой код	Размер (I.D./O.D.)
Бирюзовый	6 (1,1мм / 2,0мм)
Голубой	8 (1,7мм / 2,7мм)
Черный	10 (2,3мм / 3,3мм)
Белый	12 (2,8мм / 4,0мм)
Зеленый	14 (3,3мм / 4,7мм)
Оранжевый	16 (3,8мм / 5,3мм)
Красный	18 (4,5мм / 6,0мм)
Желтый	20 (5,1мм / 6,7мм)
Синий	22 (5,6мм / 7,3мм)
Фиолетовый	24 (6,2мм / 8,0мм)
Розовый	26 (6,9мм / 8,7мм)

Marking and labeling are made in Ukrainian and comply with all requirements of standards.

The product is registered in Ukraine and allowed to use in medical practice.

A doctor should fix  
his mind on curing



IGAR Ltd. [www.igar.ua](http://www.igar.ua)  
20, Metalistiv str., Kiev, 03057  
phone/fax: 044-585-34-55