

Gypsum bandage «IGAR»

Certificate of State Registration

№9817/2010 issued 18.10.2010



Technical characteristics:

- Excellent permeability and plasticity due to porous surface. Less loss of gypsum. Good water-proof characteristics after solidification.
- Density of supporting gauze not less than 24 g/m², of plaster - not less than 340 g/m².
- Setting time from 4 to 8 minutes.
- Time of complete solidification from 24 to 48 hours depending on size and thickness of plaster cast.

Sizes in our assortment: 10x270cm, 15x270cm, 20x27cm.

Usage: as fixing material in surgery and orthopaedy, for fixation of all type of bone fractures, orthopaedic plastic, first-aid fixation, trauma of soft tissues and joint disases.

Marking and labeling are made in Ukrainian and comply with all requirements of standards.

**A doctor should fix
his mind on curing**



**IGAR Ltd. www.igar.ua
20, Metalistiv str., Kiev, 03057
phone/fax: 044-585-34-55**
